

Open Playback Recorder

SUMMARY

Product Information

Product Name

Open Playback Recorder

Device Category

- | | |
|--|--|
| <input type="checkbox"/> Adapted Toys | <input type="checkbox"/> Kits |
| <input type="checkbox"/> Aids for Daily Living (ADL) | <input type="checkbox"/> LipSyncs |
| <input type="checkbox"/> Assistive Switches | <input type="checkbox"/> Mounting |
| <input checked="" type="checkbox"/> Communication Aids (AAC) | <input type="checkbox"/> Recreation and Leisure |
| <input type="checkbox"/> Computer Access | <input type="checkbox"/> Seating and Positioning |
| <input type="checkbox"/> Environmental Controls | <input type="checkbox"/> Switch Interfaces |
| <input type="checkbox"/> Gaming | <input type="checkbox"/> Writing Aids |
| <input type="checkbox"/> Keyguard | |

User Value Statement

The Open Playback Recorder is an open-source device that can record three lists of voice messages that its user can playback via the trigger of an accessible button connected through a 3.5 mm mono jack. It aims to assist users with communication difficulties by giving them alternative methods to engage in conversation. This device has functionality similar to AbleNet's Big Mack or Step by Step.

Designer

Design by Neil Squire Society/Makers Making Change

Device Info

Overview

This device can record and store three separate lists of messages, which can then be sequentially played back using a 3.5mm jack or an onboard button. It is designed for users with difficulty communicating, and can allow them to choose from a list of messages for three separate situations.

Disability Type

Select one or more disability types:

- | | |
|--|------------------------------------|
| <input type="checkbox"/> Agility / Dexterity | <input type="checkbox"/> Arthritis |
|--|------------------------------------|



© 2022 by Neil Squire Society

This work is licensed under the CC BY SA 4.0 License: <http://creativecommons.org/licenses/by-sa/4.0>

Open Playback Recorder

SUMMARY

- | | |
|------------------------------------|---|
| <input type="checkbox"/> Cognitive | <input checked="" type="checkbox"/> Other |
| <input type="checkbox"/> Hearing | <input type="checkbox"/> Pain |
| <input type="checkbox"/> Mobility | <input type="checkbox"/> SCI |
| <input type="checkbox"/> Mobility | <input type="checkbox"/> Vision |

Disability Type Description

This device is designed for users that have difficulty with communication and need a communication aid. This device can store three different sets of messages for use in three different situations.

How To Use

Setup

Connect two assistive switches, one to the jack labeled PLAY, and one to the jack labeled LEVEL.

Power On

To power on the device, flip the switch on the top of the device from OFF to ON.

Message Recording

To record one or more messages, the device needs to be in Record Mode. The device will record messages to the currently selected level.

Step 1. To engage “Record Mode” hold down the REC button for 2 seconds. When Record Mode is engaged, the REC LED will turn red and remain solid.

Note: When Record Mode is engaged, all previous messages on that level will be erased.

Step 2. To start message recording, press and hold the PLAY button or assistive switch connected to the play jack. The REC LED will flash red. Record the message into the microphone and release the button to end and save the message. The RED LED will once again turn red and remain solid.

Step 3. To record additional messages, repeat Step 2.

Note: These messages will later be played back in the same order they were recorded.

Step 4. When the desired messages have been recorded, press the REC button to exit Record Mode. The red REC LED will turn off.

Message Playback

To play back a message, press an assistive switch connected to the play jack. This will play a message from the current recording level. The play button on the device itself will also allow a secondary user to play a message. Continuing to press the button will cycle through all the messages on that level.



Open Playback Recorder

SUMMARY

To change the message level, press an assistive switch connected to the level jack. The level shift button on the device will also allow a secondary user to change the message level. The current level will be indicated by the three blue LEDs on the left side of the device.

Estimated Cost

The estimated material cost of the device:

- ☐ \$0 - \$10
- ☐ \$11 - \$25
- ☐ \$26 - \$50
- ☒ \$51 - \$100
- ☐ \$101 - \$250
- ☐ \$250+

Attribution

Design by Neil Squire Society/Makers Making Change

Maker Info

Project Skills

- | | |
|---|---|
| <input checked="" type="checkbox"/> 3D Printing | <input type="checkbox"/> Other |
| <input type="checkbox"/> Custom PCB | <input checked="" type="checkbox"/> Software |
| <input checked="" type="checkbox"/> Electronics | <input checked="" type="checkbox"/> Soldering |
| <input type="checkbox"/> Laser Cutting | <input type="checkbox"/> Woodworking |
| <input type="checkbox"/> Mechanics | |

Skills Description

The primary skill used in this build is soldering. The circuit is constructed on a protoboard, and all the components need to be connected by wires.

Tools Needed

- | | |
|---|--|
| <input checked="" type="checkbox"/> 3D Printer | <input type="checkbox"/> Laser Cutter |
| <input checked="" type="checkbox"/> Common Hand Tools | <input checked="" type="checkbox"/> Soldering Iron |
| <input type="checkbox"/> Common Power Tools | <input type="checkbox"/> Specialized Tooling |



Open Playback Recorder

SUMMARY

Print time (hrs)

24 hrs

Assembly time (hrs)

5 hrs

Build Instructions

The main work on this build involves soldering together two circuits on protoboards, one for each half of the enclosure, and connecting them together with Dupont wires.

Download Link

<https://github.com/makersmakingchange/Open-Playback-Recorder/archive/refs/heads/main.zip>

Project Link

<https://github.com/makersmakingchange/Open-Playback-Recorder>

License

License

- Everything needed or used to design, make, test, or prepare the Open Playback Recorder is licensed under the [CERN 2.0 Weakly Reciprocal license \(CERN-OHL-W v2\)](#) or later .
- All software is under the [GNU General Public License v3.0 \(GPL-3.0\)](#).
- Accompanying material such as instruction manuals, videos, and other copyrightable works that are useful but not necessary to design, make, test, or prepare the Open Playback Recorder are published under a [Creative Commons Attribution-ShareAlike 4.0 license \(CC BY-SA 4.0\)](#).